

midas_to_phiat Documentation

Sam Porter (wporter@triumf.ca)

University of British Columbia

TRIUMF, on behalf of the TITAN Collaboration

Overview

- **Language:** Python 3
- midas_to_phiat.py converts MIDAS files output by the MPET DAQ to readable .csv files for input into PhiAT.m
- **Input:** A directory address (i.e. address to a folder) of MIDAS Files (ex: /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_CAENOffset)
- **Output:**
 - .csv files containing position, time-of-flight and TDC trigger information
 - .csv files containing approximate count rates
 - One .csv file containing an address list of converted MIDAS files, their accumulations times (Taccs) and reference files assignments

BEFORE YOU CONTINUE! -- titan_data

- midas_to_phiat.py requires you download the **titan_data** package
 - This can be downloaded here: https://bitbucket.org/ttriumfdaq/titan_data/src/master/
- On Line 20: Change the address in `sys.path.append` from `'/Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/titan_data'` to the local address of titan_data on your computer (i.e. `'/local/address/to/titan_data'`)

```
sys.path.append('/Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/titan_data') # Address of titan_data module. NEW USERS NEED TO CHANGE THIS! #
```

A Few Other Initial Settings...

- **all_one_file:** If TRUE, resulting output data .csv files correspond to each MIDAS file
 - If FALSE, output .csv files correspond to each **loop** within in each MIDAS file (see [TITAN DAQ documentation](#))
 - **Default is TRUE** for PI-ICR analysis
- **CAEN** and **VT2:** Set TRUE whichever TDC is being used to read out data; this is indicated in the MPET DAQ under Scan
 - **Default is CAEN = TRUE**
 - Note: Cannot both be true for data converting in script -- this will trigger an error

| Misc settings | |
|---|-------------------------------------|
| Enable VT2 readout | <input checked="" type="checkbox"/> |
| Enable LRS1190 readout | <input checked="" type="checkbox"/> |
| Enable Caen V1290 readout (More settings) | <input checked="" type="checkbox"/> |

```
all_one_file = True # If TRUE, entire .mid becomes one .csv. If FALSE, e
i.e. scanned over voltages between events, change to False)
CAEN = True # If True, output data files are only CAEN TDC data
VT2 = False # If True, output data files are only VT2 TDC + LRS data
```

A Few Other Initial Settings...for Test Data

- **testing:** If you are using test data from MIDAS that **DOES NOT** include a reference file (i.e. a Tacc of 0), set to TRUE
 - Automatically assigns the last file a Tacc of 0, as needed to run through the script
- Otherwise, set to FALSE
 - This is the default for normal data

Getting Started

- midas_to_phiat.py is automatically run as part of PhIAT.m (see PhIAT.m documentation)
 - However, it can be run on its own as well:

```
# NOTE: To test run this script alone: python3 midas_to_phiat.py /path/to/directory/of/files
```

- The **input** is an address to a directory of MIDAS files (.mid or .mid.lz4)
 - The files in this directory are gathered together via **glob**

Data

- For every event (i.e. ion hit on the detector), four pieces of data are relevant:
 - **X Position** [mm]: X coordinate of ion hit on detector (i.e. delay line anode)
 - **Y Position** [mm]: Y coordinate of ion hit on detector (i.e. delay line anode)
 - **Time-of-Flight** [s]: Time-of-flight between ejection from the trap to ion hit on MCP
 - **Trigger** [int]: Number of trigger an ion hit came during (starting from 0 at file's beginning) -
- CAEN 25ps ONLY
 - **Timestamp** [s]: Timestamp of event relative to initial trigger time -- VT2 ONLY

Data -- CAEN 25ps TDC

- **X/Y Position:** Determined via timing pulses in delay line anode
 - Correspond to `caen_tdc_parsed.pos_x_mm` (or `pos_y_mm`) in `../titan_data/mpet/__init__.py`
- **ToF:** Time between trap ejection and ion contact with MCP
 - Corresponds to `caen_tdc_parsed.mcp_tof_secs`
- **Trigger:** Trigger number at which an ion event occurred
 - Corresponds to `caen_tdc_parsed.trigger_count`

Data -- VT2 TDC

- **X/Y Position:** Determined via LeCroy 1190 Position Analyzer
 - Correspond to `pos_data.x` (or `y`) in `../titan_data/mpet/__init__.py`
 - **NOTE: The LC1190 DOES NOT give reliable position data...**
- **ToF:** Time between trap ejection and ion contact with MCP
 - Corresponds to `vt2_data.tof_secs`
 - **NOTE:** Since this is raw data (i.e. each ion hit does not yet correspond to one event), need to specify a `channel_ids` to get just one trigger value per ion hit
- **Timestamp:** Timestamp of event relative to initial trigger time
 - Corresponds to `vt2_data.time_secs`
 - **NOTE:** Since this is raw data (i.e. each ion hit does not yet correspond to one event), need to specify a `channel_ids` to get just one trigger value per ion hit

```
if tdc.channel_ids == [0]:  
    event_tof_data_vt2.append(tdc.tof_secs) # ToF from trap ejection to initial MCP contact (1st ToF)  
    event_timestamp_data_vt2.append(tdc.time_secs) # Timestamp corresponding to initial MCP contact
```

Ion-Ion Interaction Data Cuts

- We want to minimize any ion-ion interaction effects in the trap on results, so we discard any data where $> N$ total ions were in the trap together
 - N corresponds to `ion_ion_cut_min`
- In practice, we throw away data with $> N$ events in the **same trigger window**
- **ONLY** works with CAEN 25ps TDC Data

Determining Ion Rates

- We want to determine the rate at which ion events are occurring
- We determine a gating time t (**gating_rate** [ms]) and determine the number of events that occur within each gating time window across the file
 - i.e. events between $0-t$, events between $t-2t$, etc.
- We approximate the *global event time* for each event as the total trap time (**trap_time**) times the trigger number of that event
- We record the *approximate global gate time*, the *number of events*, and the *rate* of those events ($\text{num_events}/t$) for each gating time window

Ion Rate .csv

- A .csv file is printed with the recorded ion rate data:
 - Column 1: Approximate Global Gate Time [ms]
 - (gating window number)* t
 - Column 2: Number of Events in each Window [int]
 - Column 3: Rate of Incoming Events in each Window [events/ms]
 - $\text{num_events}/t$
- Saved as **original_MIDAS_filename.mid_ion_rate.csv**

| | A | B | C |
|----|------|---|------|
| 1 | 50 | 0 | 0 |
| 2 | 100 | 0 | 0 |
| 3 | 150 | 0 | 0 |
| 4 | 200 | 2 | 0.04 |
| 5 | 250 | 0 | 0 |
| 6 | 300 | 0 | 0 |
| 7 | 350 | 2 | 0.04 |
| 8 | 400 | 0 | 0 |
| 9 | 450 | 0 | 0 |
| 10 | 500 | 0 | 0 |
| 11 | 550 | 2 | 0.04 |
| 12 | 600 | 0 | 0 |
| 13 | 650 | 0 | 0 |
| 14 | 700 | 2 | 0.04 |
| 15 | 750 | 0 | 0 |
| 16 | 800 | 0 | 0 |
| 17 | 850 | 0 | 0 |
| 18 | 900 | 2 | 0.04 |
| 19 | 950 | 0 | 0 |
| 20 | 1000 | 0 | 0 |

Main Data .csv

- A .csv is printed with the data described on slides 7-9
 - Column 1: X Position [mm]
 - Column 2: Y Position [mm]
 - Column 3: Time-of-Flight [s]
 - Column 4: Trigger Number [int] (CAEN) OR Timestamp [s] (VT2)
- Saved as **original_MIDAS_filename.mid_.csv**
- Main Data and Ion Rate CSVs are created for each MIDAS file in the input directory

| | A | B | C | D |
|----|------------|------------|----------|----|
| 1 | 1.62760417 | 1.62760417 | 2.44E-05 | 0 |
| 2 | 3.25520833 | 3.25520833 | 4.65E-05 | 0 |
| 3 | 1.62760417 | 1.62760417 | 2.44E-05 | 1 |
| 4 | 3.25520833 | 3.25520833 | 4.65E-05 | 1 |
| 5 | 1.62760417 | 1.62760417 | 2.44E-05 | 2 |
| 6 | 3.25520833 | 3.25520833 | 4.65E-05 | 2 |
| 7 | 1.62760417 | 1.62760417 | 2.44E-05 | 3 |
| 8 | 3.25520833 | 3.25520833 | 4.65E-05 | 3 |
| 9 | 1.62760417 | 1.62760417 | 2.44E-05 | 4 |
| 10 | 3.25520833 | 3.25520833 | 4.65E-05 | 4 |
| 11 | 1.62760417 | 1.62760417 | 2.44E-05 | 5 |
| 12 | 3.25520833 | 3.25520833 | 4.65E-05 | 5 |
| 13 | 1.62760417 | 1.62760417 | 2.44E-05 | 6 |
| 14 | 3.25520833 | 3.25520833 | 4.65E-05 | 6 |
| 15 | 1.62760417 | 1.62760417 | 2.44E-05 | 7 |
| 16 | 3.25520833 | 3.25520833 | 4.65E-05 | 7 |
| 17 | 1.62760417 | 1.62760417 | 2.44E-05 | 8 |
| 18 | 3.25520833 | 3.25520833 | 4.65E-05 | 8 |
| 19 | 1.62760417 | 1.62760417 | 2.44E-05 | 9 |
| 20 | 3.25520833 | 3.25520833 | 4.65E-05 | 9 |
| 21 | 1.62760417 | 1.62760417 | 2.44E-05 | 10 |
| 22 | 3.25520833 | 3.25520833 | 4.65E-05 | 10 |
| 23 | 1.62760417 | 1.62760417 | 2.44E-05 | 11 |
| 24 | 3.25520833 | 3.25520833 | 4.65E-05 | 11 |

Pairing Final and Reference Files

- For each final file ($T_{acc} > 0$), we need to pair it with the reference file ($T_{acc} = 0$) *closest to it in time*
- We determine the midpoint time of each file (i.e. $\text{file_start} + (\text{file_end} - \text{file_start})/2$)
 - file_start is **event_time[0]** and file_end is **event_time[-1]** from titan_data/mpet/__init__.py
- We find the difference between the midpoint time of the final file and all reference files, and pair the actual file with the reference file of *smallest* difference

```
for ref_time in ref_time_list:
    diff = np.abs(ref_time - time)
    diff_list.append(diff)

min_index = diff_list.index(min(diff_list)) # Finding the number of ref file closest in time to each actual file
minindex_list.append(min_index+1)
```

File List .csv

- A .csv is printed with:
 - Column 1: Addresses of Main File CSVs
 - Column 2: Accumulation Times of each File [s]
 - Column 3: Reference File Assignment for each Final File
 - If file is a reference file, this is blank
 - Numbering of associated references correspondings to ordering of references in first column
- Saved as **inputDirectoryName_firstTaccFile_List.csv**
 - i.e. Test_CAENOffset_0.03File_List.csv
 - Saved OUTSIDE the input directory

File List .csv Example

| | A | B | C |
|---|--|------|---|
| 1 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run00103.mid_.csv | 0.03 | 1 |
| 2 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001032.mid_.csv | 0.03 | 1 |
| 3 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001033.mid_.csv | 0.03 | 1 |
| 4 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001034.mid_.csv | 0.03 | 1 |
| 5 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001035.mid_.csv | 0.03 | 1 |
| 6 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001036.mid_.csv | 0.03 | 1 |
| 7 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001037.mid_.csv | 0.03 | 1 |
| 8 | /Users/wsporter/Documents/Physics_Research/TITAN/PIICR_Analysis/Testing/Test_Set1/run001038.mid_.csv | 0 | |
| 9 | | | |

- All reference files will come after all final files

- More on **titan_data** and the primary script referenced in this code (`../mpet/__init__.py`) can be found:
 - https://bitbucket.org/ttriumfdaq/titan_data/src/master/
 - [TITAN DAQ documentation](#)
- Please contact Sam Porter (wporter@triumf.ca) if you have any questions